

# PACE Model 300

THREE CHAMBER TEMPORARY  
CARDIAC PACEMAKER



The World's first external triple-chamber pacemaker for temporary bi-ventricular (CRT) pacing

 **OSYPKA**  
M E D I C A L

# Temporary Stimulation

Please refer to  
User's Manual for  
indications of use!

## PACE Model 300 Highlights

- Three dedicated pacing terminals
- Resynchronisation in patients with low ejection fraction
- Constant voltage output for reliable cardiac capture (similar to permanent pacemakers)
- Variable setting of AV- and VV-Delay
- Auto Sense function automatically tracks P/R wave amplitudes and adjust atrial and ventricular sensitivities accordingly
- P and R Wave amplitudes measured upon pressing PAUSE key
- Interface with Pacemaker Clinic™ for automatic CRT optimization
- Quick and easy change of battery: no need to worry about battery orientation (one single 9 V battery)
- Lifetime service



P/N 2005

## Specifications

## PACE Model 300

Modes	DDD (BV, R, L), VVI (BV, R, L), AAI, VDD (BV, R, L) D00 (BV, R, L), V00 (BV, R, L), A00, DVI (BV, R, L), DAI (BV, R, L), VAT (BV, R, L), AAT, DDD+AT, DAT (BV, R, L), DDI (BV, R, L)	
Basic Pacing Rates	30 – 220 ppm	
Rapid Atrial Pacing Rates	70 – 1,000 ppm	
Output Amplitude	0.1 – 18 V	
Pulse Width	0.05 – 1.5 ms	
V-V Delay	5 – 80 ms	
Sensitivity Threshold	Atrium : 0.2 – 20 mV	Ventricle: 1.0 – 20 mV
Refractory Period	250 ms	
Dimensions	21.7 cm x 9.6 cm x 5.1 cm (8.9" x 3.8" x 2.0")	
Weight incl. Battery	490 g (17.3 ounces)	
Battery Type	9V Type 6LP3146 Alkaline	
Battery Life	Alkaline: Up to 5 days (72 ppm, 8.0 V)	



Osypka Medical GmbH

Albert-Einstein-Strasse 3  
12489 Berlin, Germany  
Phone: +49 (30) 6392 8300  
Fax: +49 (30) 6392 8301  
E-Mail: mail@osypkamed.com  
www.osypkamed.com

United States of America:  
Osypka Medical, Inc.

dba. Cardiotronic  
La Jolla, CA 92037, USA  
Phone: +1 (858) 454 0600  
Fax: +1 (858) 454 0640  
E-Mail: mail@cardiotronic.net  
www.cardiotronic.net

Copyright © 2017 Osypka Medical

